



The
Body Images
Mustang



*Penis
Enlargement
System*

Instructions for Use

INSTRUCTIONS FOR USE
for
THE MUSTANG PENIS ENLARGING SYSTEM

1. Trim or shave the pubic hair at the base of the penis.
2. A warm bath, shower or Jacuzzi before a pump session is recommended
3. Massage and stretch the penis and scrotum 3 to 5 minutes prior to a pump session.
4. Lubricate generously, the penis head, shaft and base. For a secure seal, apply a lubricant directly on the pubic hairs. Always use a good quality lubricant (i.e Body Images Extra Sensual Love Lube). Re-massage and re-lubricate the penis and scrotum at least every 15 to 20 minutes.
5. Do not ejaculate several hours prior to or during a pump session. This allows maximum expansion without swelling the penis.
6. Be semi to near full erection while using vacuum pressure and try to be sexually aroused throughout the pump session.
7. Place the cylinder over the penis with the valve at the top shut and the tubing connected to the vacuum pump.
8. Use minimum pressure only, you should never feel pain or discomfort. If swelling, discoloration or pain occur reduce the pressure immediately and discontinue pump usage. Do not try to fill the cylinder. When the penis becomes tight around the base of the cylinder and/or if it is difficult to remove, **a larger cylinder is required. The purposes of the pressure is to create an erection slightly larger than your normal erection, and NOT TO SWELL YOUR PENIS.**
9. Pump sessions should last 20 to 30 minutes during the first week. Increase the sessions 5 to 10 minutes each week, until you reach a pumping time of one hour per session.
10. When using two cylinders, alternate each session with different cylinder sizes for 6 days (on day 1, 3 and 5, use the small cylinder . Day 2, 4 and 6 use the large cylinder, then rest on days 7 and 8). Never use the same cylinder twice within 24 hours. This allows time for the tissue to relax and rebuild.
11. Always use clean cylinders. Wash with soap and water – do not share cylinders with friends.
12. Wash your penis and urinate after each pump session. Itching may occur if you do not clean all lubricated areas thoroughly.

The vacuum can be released by pressing the button on the side of the pump or by unscrewing the valve on the top of the cylinder slightly.

‘THE TWO CYLINDER METHOD’

Two expansion cylinders are required for optimum enlargement. One cylinder covers the penis, the second cylinder covers the penis and scrotum.

Both Cylinders are required for maximum penis enlargement. The first or penis cylinder will thicken as well as lengthen the penis. The penis scrotum cylinder or second cylinder pulls at the suspensory ligaments and gives additional length.

The second cylinder is used to increase the penis length. It is not a scrotum cylinder. Although the scrotum can become fuller in appearance, the testicle cannot get bigger using the vacuum pump.

Warning: The vacuum pump device can maintain high levels of pressure which may cause injury. Use minimal pressure only and read the instructions fully. The device is not intended for medical use. Individuals with a medical or physical condition or any persons with implants should not use the vacuum pump and should consult a doctor prior to use. Usage is done at your own risk. This is a recreational device sold only as a novelty item.